

MR Uyumlu Aspiratör ve Aksesuarları

MRI Compatible Suction Unit and Accessories

- 50 PSI ayarlı diyafram ve 0-15 LPM flowmetre.
- Herhangi bir 50 PSI oksijen kaynağıyla kullanılabilir.
- Kullanım için elektrik kaynağı veya pil gerektirmez.
- 0-350 mm Hg çekiş ve 21 LPM akış hızı üretir.
- Yüksek darbeye dayanıklı toplama kabı
- Toplama kabı kirliliğini ortadan kaldıran hidrofobik filtreye sahiptir.
- Çift 50 PSI kontrol valf çıkışı, aynı anda çekme ve oksijen kullanımına izin verir.
- Hastanın taşınması veya acil arabalarına takmak için idealdir.
- Comes with 50 PSI preset diaphragm regulator & 0-15 LPM flowmeter.
- Runs off any 50 PSI oxygen source.
- No electricity or batteries needed to operate.
- Produces 0-350 mm Hg of suction & flowrate of 21 LPM.
- High impact resistant collection canister.
- Collection canister has built in hydrophobic filter which eliminates contamination.
- Dual 50 PSI check valved outlets allow use of suction and oxygen simultaneously.
- Great for patient transport or attaching onto crash carts.

